# Instructions for Use

380157-A

English

## LZ-PEPSINOGEN II

REF V-IZ53

### INTENDED USE

For measurement of Pepsinogen II in serum or plasma.

### INTRODUCTION

Pepsinogen is an inactive zymogen of pepsin, a proteolytic enzyme found in gastric secretions. It is immunochemically classified into two types, pepsinogen I (PG I), and pepsinogen II (PG II). PG I is secreted by the chief cells and accessory cells in the gastric fundus and can also be detected in the blood and urine. Besides the gastric fundus, PG II is found in the glands of cardia, pylorus, and duodenum, and can also be detected in the blood. Since the blood levels of PG I and II reflect primarily the number of chief cells and accessory cells in the gastric mucosa, it is thought that assays of PG I and PG II in the blood, as well as the I/II ratio, may be useful in the evaluation of the degree of atrophic gastritis, diagnosis of peptic ulcer disease, and assessment of gastric acidity.

Chronic atrophic gastritis is now thought to be a premalignant lesion. It is known that lesions with a high degree of atrophy are associated with a higher incidence of stomach cancer, and it is also likely that blood levels of PG I and PG II, which reflect the degree of severity of gastric mucosa atrophy, are useful as a screening tool in populations at high risk for the development of stomach cancer.

LZ-PEPSINOGEN II is a reagent developed for highly sensitive and accurate measurement of PG II. This measurement utilizes a latex agglutination reaction, and the change in turbidity caused by this reaction is measured optically to determine the concentration.

### PRINCIPLE OF THE METHOD

This method is an optical measurement method using a latex agglutination reaction and an automated analyzer.

The latex reagent is prepared by binding anti-PG II antibodies to the surface of the latex particles. When this reagent is mixed in a cell to react with the sample, the anti-PG II antibodies that were bounded to the latex particles react with the PG II in the sample, causing agglutination.

This reaction is then measured as a change in the turbidity, with the amount of the change increasing in proportion to higher concentrations of PG II in the sample.

Measurement using LZ-PEPSINOGEN II applies this principle to find a calibration curve from standards of known concentration. The amount of PG II in the sample is then found relative to this standard.

### CONTENTS OF THE KIT

1. Reagent-1 .....60mL, 1 vial  
(Contains 50mmol/L of N-2-Hydroxyethylpiperazine-N'-2-ethanesulfonic acid (HEPES))
2. Reagent-2 ..... 20mL, 1 vial  
(Contains 50 vol% of latex sensitized with anti Pepsinogen II mouse monoclonal antibodies)

### WARNINGS AND PRECAUTIONS

1. For in vitro diagnostic use only.
2. Use the fresh serum or plasma. When samples are stored, they should be kept at -20°C. Repeated freezing and thawing of sample should be avoided.
3. Store the reagents under the designated conditions. Do not use reagents that have passed their expiration date.
4. Mix the latex reagent before using by gently inverting the vial several times.
5. Measurement errors may result if bubbles are present on the surface of the sample after it is dispensed into the sample cup. Therefore remove all bubbles. If fibrin is present in the sample cup, remove it. Fibrin can cause clogging of the sample nozzle.
6. Create a calibration curve for each day of measurement. Also be sure to create a new calibration curve when a reagent from a different vial or lot is used.
7. The test sample may be contaminated with the HB virus, HCV, HIV, or other pathogens. Therefore use caution when handling.
8. If the sample concentration exceeds the measurement range, dilute with a normal saline solution or similar solution and perform measurement again.
9. If a reagent enters the eye or mouth, rinse it out with large volumes of running water, and perform other required first aid. If necessary, seek medical attention.
10. Use the reagents as quickly as possible after they are opened. If they are to be stored, be sure to close the caps and store them using the prescribed method.
11. There is the danger of infection from all tools, reagents, and reagent containers that contact the sample. Disinfect them using an autoclave or other means, or soak them in hypochlorous acid or other disinfectant solution.  
Example of treatment: Soak for 60 minutes or longer in a sodium hypochlorite solution (available chlorine concentration 1000 ppm or greater). (Neutralize any substances that contain acids before soaking.) Alternatively, treat in an autoclave at 121°C for 20 minutes. (Do not treat in this way any items to which sodium hypochlorite has adhered.)
12. Dispose of used reagents and containers as medical waste in accordance with applicable regulations.
13. If the product is used in any way other than that specified here, the reliability of measurement results cannot be guaranteed. Be sure to follow the procedure.
14. A clinical diagnosis based on the measurement results must be a comprehensive judgment made by the attending physician, including factors such as clinical symptoms and other test results.

### SAMPLE COLLECTION

1. Use serum or plasma as test sample (specimen).
2. Collect samples in usual manner and test them while fresh.
3. If test samples are to be stored over a long period of time, preserve them frozen at -20°C or lower (avoid repeated freezing and thawing).
4. If a frozen stored test sample is to be used, thaw it at room temperature and then mix it by inversions before running the test.
5. Dialyzed blood-derived serum samples, in which fibrin is likely to be separated, should be defibrinated before used in the test.

### TEST PROCEDURE

#### 1. Preparation of reagents

- 1) Reagent-1: Use the Reagent-1 as-is.
- 2) Reagent-2: Use the Reagent-2 as-is.
- 3) Standards: Separately obtain the LZ-PG II Standard (REF V-IZ96) as the standards.

#### 2. Measurement procedure

- 1) Follow the instructions in the instruction manual to operate the automated analyzer.
- 2) Dispense the sample and standards into sample cups, and set them in the designated positions.
- 3) Set Reagent-1 and Reagent-2 in the designated positions.
- 4) Enter the parameters into the analyzer.
- 5) Press the analyzer "Start" key to begin analysis and output the measurement results (by printout or other means).

### EXPECTED REFERENCE VALUES

When pepsinogen I and II are measured simultaneously,

- Positive :Pepsinogen I is below 70ng/mL  
and the ratio of Pepsinogen I / II is below 3
- Strong Positive :Pepsinogen I is below 30ng/mL  
and the ratio of Pepsinogen I / II is below 2

### INTERFERING SUBSTANCES

Almost no effect on the measurement value was found from conjugated bilirubin (20mg/dL), free bilirubin (20mg/dL), hemoglobin (500mg/dL), chyle (2240 formazine turbidity units), and rheumatoid factor (RF positive sample 500 IU/mL). As an anticoagulant almost no effect on the measurement value was found from EDTA · 2Na (200mg/dL), sodium citrate (500mg/dL) and heparin sodium (50mg/dL).

### INTERNAL QUALITY CONTROL

Each laboratory should establish a quality control program to monitor the performance of LZ-PEPSINOGEN II. It is recommended to use the following materials for the quality control in your laboratory.

- PG-Control Low (REF V-IZ97)  
PG-Control High (REF V-IZ98)

### PERFORMANCE CHARACTERISTICS

#### 1. Sensitivity

When standard solution of 2ng/mL Pepsinogen II was measured in comparison to the blind reagent test, the change in absorbance that was obtained was 0.002 or more.

#### 2. Accuracy

When control samples of known concentration were measured, the value obtained was within ±15% of indicated value.

#### 3. Within-run reproducibility

When the same sample was measured 10 times within the same run, the coefficient of variation (CV) for the values obtained was 10% or less.

#### 4. Measurement range

1 – 100ng/mL

### PRODUCT CODE, PRODUCT NAME & STORAGE

Product code	Product name	Contents	Storage
V-IZ53	LZ-PEPSINOGEN II	1 × 60mL 1 × 20mL	2-8°C
V-IZ96	LZ-PG II Standard	6 × 1mL	2-8°C
V-IZ97	PG-Control Low	10 × 1mL	2-8°C
V-IZ98	PG-Control High	10 × 1mL	2-8°C

### REFERENCE

1. P. Correa et al.: Cancer Survey, 2: 437-450, 1983.

WELLKANG TECH CONSULTING  
Suite B, 29 Harley Street, London W1G 9QR, UK

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EIKEN CHEMICAL CO., LTD.