# Instructions for Use

380158-B

English

## LZ-PG I Standard

REF V-IZ95

### INTENDED USE

The product is intended for use as the standards when using automated analyzers to measure Pepsinogen I in serum or plasma.

### REAGENTS

Reagents are stable until the date printed on the label assuming the container remains unopened at a storage temperature of 2-8°C.

LZ-PG I Standard                      6 × 1 mL, Ready to use.

### WARNINGS AND PRECAUTIONS

1. For *in-vitro* diagnostic use only.
2. Do not freeze product.
3. Do not use expired product.
4. Air bubbles on the surface of the product after being poured into sample cups can result in error measurement, and should be removed.
5. LZ-PG I Standard is manufactured from human-derived components. Only components testing negative for HBs antigens, HIV (HIV-1 and HIV-2) antibodies and HCV antibodies is used, but the presence of pathogens and infection risk cannot be completely ruled out. Therefore, the product should be handled in the same manner as patient samples to ensure safety.
6. If a reagent enters the eye or mouth, rinse it out with large volumes of running water, and perform other required first aid. If necessary, seek medical attention.
7. Use the reagents as quickly as possible after they are opened. If they are to be stored, be sure to close the caps and store them using the prescribed method.
8. Dispose of used reagents and containers as medical waste in accordance with applicable regulations.

### PROCEDURE

Create the calibration curve in the manner described by automatic analyzer manual.

### PRODUCT CODE, PRODUCT NAME & STORAGE

Product code	Product name	Contents	Storage
V-IZ95	LZ-PG I Standard	6 × 1mL	2-8°C

### CONCENTRATION

Standard Number	Concentration (ng/mL)
1	0.0
2	10.0
3	25.0
4	70.0
5	130.0
6	200.0

**Advena Ltd.**  
Tower Business Centre, 2nd Flr., Tower Street,  
Swatar, BKR 4013 Malta

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EC REP

Advena Ltd.

Tower Business Centre, 2nd Flr., Tower Street,  
Swatar, BKR 4013 Malta

EIKEN CHEMICAL CO., LTD.

4-19-9, Taito, Taito-ku, Tokyo 110-8408, Japan