Instructions for Use

380160-B

English

PG-Control Low

REF V-IZ97

PG-Control High

REF V-IZ98

INTENDED USE

The product is intended for use as an internal quality control fluid in medical laboratories making use of automated analyzers to measure Pepsinogen I and Pepsinogen II.

REAGENTS

Reagents are stable until the date printed on the label assuming the container remains unopened at a storage temperature of 2-8°C.

WARNINGS AND PRECAUTIONS

1. For *in-vitro* diagnostic use only.
2. Do not freeze product.
3. Do not use expired product.
4. Air bubbles on the surface of the product after being poured into sample cups can result in error measurement, and should be removed.
5. PG-Control is manufactured from human-derived components. Only components testing negative for HBs antigens, HIV (HIV-1 and HIV-2) antibodies, and HCV antibodies is used, but the presence of pathogens and infection risk cannot be completely ruled out. Therefore, the product should be handled in the same manner as patient samples to ensure safety.
6. Pepsinogen I and pepsinogen II concentration is displayed on the attached data sheet.
7. The product should be used with LZ-PG I Standard (REF V-IZ95) or LZ-PG II Standard (REF V-IZ96).
8. PG-Control should be stored at 2-8°C after reconstitution and used within 7 days.
9. Dispose of used reagents and containers as medical waste in accordance with applicable regulations.

PROCEDURE

1. Gently remove the rubber stopper from the vial. Add 1.0mL of deionized water.
2. Seal the vial with the rubber stopper, and mix by gently inverting the vial several times. Reconstitute the control at least 15 minutes before use.
3. Dispense an appropriate quantity of the reconstituted control into a sample cup and then operate the automatic analyzer according to its instructions.

PRODUCT CODE, PRODUCT NAME & STORAGE

| Product code | Product name | Contents | Storage |
|--------------|-----------------|----------|---------|
| V-IZ97 | PG-Control Low | 10 × 1mL | 2-8°C |
| V-IZ98 | PG-Control High | 10 × 1mL | 2-8°C |

Advena Ltd.
Tower Business Centre, 2nd Flr., Tower Street,
Swatar, BKR 4013 Malta

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Advena Ltd.

Tower Business Centre, 2nd Flr., Tower Street,
Swatar, BKR 4013 Malta

EIKEN CHEMICAL CO., LTD.

4-19-9, Taito, Taito-ku, Tokyo 110-8408, Japan

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