



Instructions for Use

380242-B

English

H.pylori Antibody Control Low

REF V-IY84

H.pylori Antibody Control High

REF V-IY85

INTENDED USE

The product is intended for use as an internal quality control fluid in medical laboratories when using automated analyzers to detect *Helicobacter pylori* antibodies.

REAGENTS

Reagents are stable until the date printed on the label assuming the container remains unopened at a storage temperature of 2-8°C.

WARNINGS AND PRECAUTIONS

1. For *in-vitro* diagnostic use only.
2. Do not freeze product.
3. Do not use expired product.
4. Air bubbles on the surface of the product after being poured into sample cups can result in error measurement, and should be removed.
5. No human-derived ingredients are contained in this product, however, the risk of infection and the possible existence of other pathogens cannot be ruled out completely. Handle it as carefully as a patient's specimen.
6. Please refer to the attached data sheet of *Helicobacter pylori* antibody concentration.
7. The product should be used with LZ *H.pylori* Antibody Calibrator (**REF** V-IY83).
8. If a reagent enters the eye or mouth, rinse it out with large volumes of running water, and perform other required first aid. If necessary, seek medical attention.
9. Use the reagents as quickly as possible after they are opened. If they are to be stored, be sure to close the caps and store them using the prescribed method.
10. Dispose of used reagents and containers should be performed as medical waste in accordance with applicable regulations.

PROCEDURE

Dispense an appropriate quantity of the control into a sample cup and then operate the automatic analyzer by following to its instructions.

PRODUCT CODE, PRODUCT NAME & STORAGE

Product code	Product name	Contents	Storage
V-IY84	H.pylori Antibody Control Low	2 × 3mL	2-8°C
V-IY85	H.pylori Antibody Control High	2 × 3mL	2-8°C



Advena Ltd.

Tower Business Centre, 2nd Flr., Tower Street, Swatar, BKR 4013 Malta



EIKEN CHEMICAL CO., LTD.

4-19-9, Taito, Taito-ku, Tokyo 110-8408, Japan

(Date of Issue: April 1, 2019)