# Instructions for Use

380340-A

English

## OC-SENSOR Sample Diluent **REF** V-PH19

### INTENDED USE

OC-SENSOR Sample Diluent is designed for use as a high concentration sample diluent for the determination of hemoglobin or calprotectin in feces with OC-SENSOR automated analyzers.

### PROCEDURE

Manually dilute high concentration samples with this product before measurement with OC-SENSOR automated analyzers [OC-SENSOR Micro, OC-SENSOR io].

Automatically perform a dilution retest in accordance with the operation manual for OC-SENSOR automated analyzers. Set a high-concentration sample, an empty sample cup and a sample cup with 2 mL (or more) of this product on a dedicated rack. Measure with OC-SENSOR automated analyzers [OC-SENSOR DIANA, OC-SENSOR PLEDIA].

### STORAGE

This product is stable up to the date printed on the label as long as it remains unopened at a storage temperature of 2-8°C.

### COMPOSITION AND DESCRIPTION

This product contains buffer solution (HEPES; N-2-hydroxyethyl piperazine-N'-2-ethanesulfonic acid). The container of this product is made of polyethylene (PE), the cap is made of polypropylene (PP), and the case is made of paper.

### WARNINGS AND PRECAUTIONS

1. For *in vitro* diagnostic use only.
2. Do not use this product for any purpose other than dilution of high concentration samples.
3. Do not freeze products.
4. Do not use expired products.
5. Air bubbles on the surface of the product after being poured into sample cups can result in error measurement, and should be removed.
6. Keep this product out of the reach of children.
7. The test sample may contain microorganisms. Therefore, please be very careful when handling them. After use, all specimens, etc. should be considered to be at risk of infection and treated accordingly.
8. Dispose of used reagents and containers as medical waste in accordance with local regulations.
9. This product contains sodium azide ( $\text{NaN}_3$  <0.1%), which is a toxic substance. Take care not to touch skin, eyes or mouth. If exposed, rinse well with water and contact your doctor for treatment advice.

### PRODUCT CODE, PRODUCT NAME & STORAGE

Product code	Product name	Contents	Storage
V-PH19	OC-SENSOR Sample Diluent	3 x 45mL	2-8°C

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Advena Ltd.  
Tower Business Centre, 2nd Flr., Tower Street,  
Swatar, BKR 4013 Malta

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**EIKEN CHEMICAL CO., LTD.**