# Instructions for Use

380273-A

English

**REF** 978000

## Loopamp™ MALARIA Pf

### Detection Kit

#### INTENDED USE

The Loopamp™ MALARIA Pf Detection Kit is a qualitative *in vitro* diagnostic test for the detection of *Plasmodium falciparum* DNA extracted from human blood samples.

#### TEST PRINCIPLES

This product is based on the nucleic acid amplification method, LAMP (Loop-mediated Isothermal Amplification), developed by Eiken Chemical Co., Ltd.

The LAMP method has the following characteristics: (1) Only one polymerase enzyme is required and the amplification reaction proceeds under isothermal conditions;<sup>1,2</sup> (2) it has an extremely high specificity because of the use of four primers recognizing six distinct regions on the target; (3) it has a high amplification efficiency and can produce a high concentration of amplified product in a short time, which makes visual or automated detection possible.<sup>3,4</sup>

The targeted DNA sequences of the *P. falciparum* (Pf)-specific primers have been designed to detect the mitochondrial DNA and have been confirmed by alignment analysis to be specific for Pf-sequences.

The test DNA solution extracted from blood samples is dispensed into a reaction tube. The strand displacement DNA polymerase, deoxynucleotide triphosphates (dATP, dCTP, dGTP and dTTP), calcein, reaction buffers and Malaria Pf-specific primers are stored in dried form in the cap of the reaction tube. These dried LAMP reagents (Malaria Pf detection reagents (dMAL Pf)) dissolve when the DNA solution is added. The reaction tube is then incubated at 65.0°C and the DNA is amplified by the strand displacement DNA polymerase in accordance with the LAMP reaction.

The detection of amplified products is based on turbidity of magnesium pyrophosphate (a white precipitate produced as a by-product of DNA amplification).<sup>3</sup> Alternatively, visual detection under UV light may be used. Before the DNA amplification, calcein contained in the reagent is in the quenched state as it is bound to manganese ions. At the start of the DNA amplification, the pyrophosphate ions that are generated out-compete the manganese ions for binding sites, and thus the calcein becomes fluorescent.<sup>4</sup>

#### CONTENTS OF THE KIT

Reagents are stable until the date on the label assuming the container remains unopened at a storage temperature of 2 – 30°C.

Malaria Pf detection reagent (dMAL Pf) ..... 2 x 48 tubes

The following reagents in dried form are contained in each reaction tube:

- Bst* DNA polymerase<sup>1</sup>
- Deoxynucleotide triphosphates
- Magnesium sulfate
- Calcein
- Manganese chloride
- Primers<sup>2</sup>

Positive control Mal (PC Mal)<sup>3</sup> ..... 1 x 1.0 mL

Negative control Mal (NC Mal) ..... 3 x 0.5 mL

30µL dropper ..... 5 x 12 droppers

\*1: *Bst* DNA polymerase derived from *Bacillus stearothermophilus* is a strand displacement DNA polymerase that lacks 5'→3' exonuclease activity.

\*2: Primers designed for the mitochondrial DNA of Pf, purified from synthesized oligonucleotides by HPLC.

\*3: PC Mal contains a product resulting from *in vitro* amplification of an artificial gene designed from the mitochondrial DNA of Pf (GenBank No.M76611).

The abbreviations of names of the following reagents, their Lot No. and

the manufacturer (EKN), are printed on the containers.

Reagents	Labelling on the tube	Code on the cap
Positive control Mal	PC Mal Lot No., EKN	PC Mal
Negative control Mal	NC Mal Lot No., EKN	NC Mal

#### WARNINGS AND PRECAUTIONS

- (1) For *in vitro* diagnostic use only.
- (2) This product is designed only for the detection of DNA of Pf parasites in blood samples of human origin. Do not use for other purposes.
- (3) When using this product, always follow this Instructions for Use.
- (4) Do not freeze the reagents.
- (5) Do not use any expired reagent.
- (6) Do not mix reagents from different Lots.
- (7) Do not replenish any reagent.
- (8) Performance of the Loopamp™ MALARIA Pf Detection Kit is dependent on operator proficiency and adherence to procedural directions. Testing should be performed by properly trained personnel strictly according to the instructions provided.
- (9) Exposure to heat, humidity and light may cause deterioration of the dMAL Pf. Remove only the required number of reaction tubes (sum of samples and controls) and re-seal the aluminium pouch immediately.
- (10) Do not remove the desiccant from the aluminium pouch. High humidity can cause deterioration of the dried LAMP reagent in the reaction tubes.
- (11) Read the instruction manual and ensure required equipment (turbidimeter or incubator) is available before commencing the procedure.
- (12) Blood samples pose a potential risk for infection. Use universal precautions to minimize biohazard.<sup>5</sup>
- (13) PC Mal and NC Mal contain a small amount of sodium azide as preservative. Sodium azide is classified as toxic. Avoid any contact with eyes, mouth, or skin.
- (14) In case of accidental contact of any reagent with eyes, mouth, or skin, immediately rinse the affected site with running water and seek medical advice.
- (15) Do not dilute or add the PC Mal to the samples. Use the PC Mal only as described in this package insert in order to avoid DNA contamination.
- (16) Store the PC Mal and any positive blood samples separately from the other kit reagents.
- (17) The cap of each reaction tube contains dMAL Pf in dried form. Do not touch the inside of the cap.
- (18) Before using the reaction tubes, check carefully to see if they have any cracks or scratches. Damaged tubes might give false results and lead to DNA contamination of the incubator and work area.
- (19) Do not expose reaction tubes to UV light before the end of the LAMP reaction. Prolonged exposure to UV light might damage the tubes and lead to false results.
- (20) When a UV light is used for visual fluorescence judgment, do not stare directly at UV light. Since UV light is harmful to the eyes, even watching for a short period may irritate eyes and cause symptoms similar to conjunctivitis. Use a glass screen or wear protective goggles/glasses/face mask whenever looking directly at the UV light.
- (21) Refer to the manual of the incubator. When the HumaLoop M or the Real-Time Turbidimeter HumaTurb C+A is used, be careful when removing the reaction tubes from the incubator to avoid burns.
- (22) Do not use PC Mal as positive control for Loopamp™ MALARIA Pan Detection Kit, Loopamp™ MALARIA Pv Detection Kit, Loopamp™ MALARIA Pan/Pf Detection Kit. Do not use PC of other kit as positive control for this kit.

#### WASTE DISPOSAL

- (1) Do not open the tubes after DNA amplification. Leave the cap closed and dispose the used tubes as medical waste by incineration or after double bagging with sealable plastic bags.
- (2) Never autoclave or reuse the reaction tubes. Amplified products will disperse and cause contamination.
- (3) The main material for the reaction tubes and reagent tubes is PP; for the reaction tube tray, PET; for the aluminium pouch, aluminium; and

for the kit case, paper.

- (4) Dispose of any other reagent, container, or lab ware in accordance with local regulations.

## SPECIMEN COLLECTION

- (1) Blood samples should be used immediately after collection.
- (2) Collect blood in a separate room from the LAMP amplification room. Aerosols containing Pf DNA can be generated during blood collection and may cause contamination.
- (3) DO NOT USE EDTA and Citrate as anticoagulant for blood collection if the result is to be read by fluorescence. The use of heparin as anticoagulant is recommended.

## MATERIALS REQUIRED BUT NOT PROVIDED

- Loopamp™ PURE DNA Extraction Kit (optional)
- HumaHeat (optional)

### For visual fluorescence detection

#### (For HumaLoop M)

- HumaLoop M

#### (For other incubator using UV light)

- Incubator (temperature accuracy:  $\pm 0.5^{\circ}\text{C}$ ; with hot bonnet)
- Heating block
- UV light or Blue LED light (wavelength: 240 to 260 nm, and 350 to 370 nm)
- Goggles/glasses or other UV-blocking eye mask

### For real-time turbidity detection

- HumaTurb C+A

### For reagent and sample mixing

- Micropipettes (10 to 100  $\mu\text{L}$ , and 20 to 200  $\mu\text{L}$ ) and pipette tips with filter
- Centrifuge for micro-tubes (optional)
- Centrifuge for eight connected tubes (optional)
- HuMax ITA (optional)

## PREPARATION OF SAMPLE DNA SOLUTION

To extract the DNA from blood samples, boil & spin and PURE method are recommended. The standard operating procedures (SOP) of both methods are provided at

<https://www.finddx.org/wp-content/uploads/2016/04/SOP-LAMP-Malaria-Aug2012.pdf>

Pay attention to the following critical points for PURE method.

- **Samples : Blood samples collected by finger prick, collected in heparin tubes or whole blood dried on filter paper**
- **Sample volume : 30 $\mu\text{L}$  (fresh blood or blood with heparin) or 6-mm blood spot punch (dried blood spot)**
- **Additive : Add 30 $\mu\text{L}$  of 334mM NaCl solution (not included in the Loopamp™ PURE DNA Extraction Kit) to the heating tube before heating**
- **Heating : For 5 minutes at 75 $^{\circ}\text{C}$**

## PREPARATION OF REAGENTS

### (1) Malaria Pf detection reagent

Remove the required number of tubes from the aluminium pouch and put them in the provided rack (sum of samples and controls).

*Note: After removing the required tubes, reseal the aluminium pouch with any unused tubes immediately.*

### (2) Negative control Mal (NC Mal)

Flick (or spin) down the tube before using, in order to collect the content at the bottom of the tube.

*Note: NC Mal should be used with every run.*

### (3) Positive control Mal (PC Mal)

Flick (or spin) down the tube before using, in order to collect the content at the bottom of the tube.

*Note: PC Mal should be used with every run.*

## MEASUREMENT PROCEDURE

### Reagent and sample mixing

- (1) Turn on the HumaLoop M or the Real-Time Turbidimeter HumaTurb C+A.
- (2) Dispense 30  $\mu\text{L}$  of extracted DNA solution into a reaction tube using the pipette and close the cap.
- (3) Dispense 30  $\mu\text{L}$  of NC Mal into a reaction tube using the pipette or dropper, and close the cap.

- (4) Dispense 30  $\mu\text{L}$  of PC Mal into a reaction tube using the pipette or dropper, and close the cap.

- (5) Flick (or spin) down all tubes to collect the solution at the bottom of the tubes.

*Note: When using the PURE device make sure the liquid level is closer to the upper line of the two lines on a reaction tube to ensure that the correct volume has been dispensed.*

- (6) Reconstitute the dried reagents in the cap by inverting the reaction tubes and collecting the DNA solution in the cap. Leave the tubes standing upside down for 2 minutes to reconstitute the dried reagents.

- (7) Invert the reaction tubes 5 times to mix the contents. Make sure that the dried reagents in the cap are fully dissolved (the solution should have a slight orange colour).

- (8) Flick (or spin) down all tubes to collect the solution at the bottom of the tubes.

## Amplification

### For visual fluorescence detection

#### (For HumaLoop M)

- (1) Check that the temperature of the HumaLoop M is 65.0 $^{\circ}\text{C}$ .
- (2) Load the reaction tubes into the HumaLoop M and press the green button to start the LAMP reaction (40 minutes at 65.0 $^{\circ}\text{C}$ ). See the HumaLoop M instruction manual for details on how to operate the incubator.
- (3) Confirm the completion of polymerase inactivation (automatically completed by the HumaLoop M). Take all reaction tubes from the HumaLoop M.

#### (For other incubator using UV light)

- (1) Set the temperature of the incubator to 65.0 $^{\circ}\text{C}$  (with hot bonnet temperature set to 10 degrees above the reaction temperature or as near to this figure as possible – temperature accuracy:  $\pm 0.5^{\circ}\text{C}$ ). Wait until the displayed temperature reaches the set value.
- (2) Load the reaction tubes and then start amplification reaction (for 40 minutes at 65.0 $^{\circ}\text{C}$ ).
- (3) 40 minutes later, inactivate the polymerase using the heating block (for 5 minutes at 80 $^{\circ}\text{C}$ , or for 2 minutes at 95 $^{\circ}\text{C}$ ) to terminate the reaction.

### For real-time turbidity detection with HumaTurb C+A (see Flow chart of the procedure)

- (1) Configure the Real-Time Turbidimeter HumaTurb C+A for detection with this product.
- (2) Check if the displayed temperature reaches 65.0 $^{\circ}\text{C}$  (allow the turbidimeter to warm up for 20 minutes before use).
- (3) Insert the reaction tubes and start measurement.
- (4) Watch the display of the turbidimeter to check the positive and negative controls for any increase in turbidity. If the turbidity increases in the positive but not in the negative control solution, amplification reaction is proceeding properly (Fig1). If this is not the case, amplification reaction may be proceeding in a wrong way. In such a case, retest affected samples from reagent preparation.
- (5) Confirm the completion of polymerase inactivation (automatically completed by the turbidimeter). Take all reaction tubes from the HumaTurb C+A and discard them without opening.

## Amplification plots by Malaria Pf detection reagent

(Analyzer : Real-Time Turbidimeter HumaTurb C+A)


Fig1 : Amplification plots for controls

Fig2 : Amplification plots for samples

## PROCEDURAL NOTES

- (1) The LAMP reaction is very sensitive and a contamination with small amounts of amplified product might lead to false positive results.
- (2) Separate the sample preparation and the amplification areas.

- (3) Clean benches with over 0.5% sodium hypochlorite before and after performing the test.
- (4) Change gloves after transferring the blood or if the gloves come in contact with the DNA solution.
- (5) When handling this product, avoid microbial contamination and nuclease contamination. Even a small amount of contamination of the reaction tube from sweat or saliva may decompose the DNA and cause a false result.
- (6) Furthermore, read the SOP<sup>6)</sup> when performing DNA extraction.
- (7) The DNA solution should ideally be used immediately after preparation.
- (8) **(For HumaLoop M or other incubator using UV light)**  
If bubbles are present, flick the tubes to release them.  
**(For Real-Time Turbidimeter HumaTurb C+A)**  
Since bubbles in reaction solution may interfere with turbidity measurement and cause a false result, avoid forming any bubble when mixing reagent and sample solution. If bubbles occur, spin or flick the tube to release them.
- (9) dMAL Pf should be fully dissolved. Any undissolved portion may influence the performance, such as causing a decrease in sensitivity.
- (10) The PC Mal contains a high copy number of control DNA. Avoid any contamination of other samples with the PC Mal. Dispense the samples and the NC Mal and close all reaction tubes before dispensing the PC Mal.
- (11) Flick (or spin) down the PC Mal tube before opening it, in order to collect the content at the bottom of the tube. Close the tube immediately after dispensing the PC Mal.
- (12) Never open the reaction tubes once the LAMP reaction has started or after completion. Be particularly careful when unloading the reaction tubes from the incubator to avoid opening the tubes accidentally.
- (13) When the HumaLoop M or the Real-Time Turbidimeter HumaTurb C+A is used, polymerase inactivation is automatically performed.
- (14) For other incubators, when visual fluorescence judgment is chosen, inactivate the polymerase (for 5 minutes at 80°C, or for 2 minutes at 95°C) before reading, or false results will be caused.
- (15) Do not reuse any amplified product in the tubes for electrophoresis or other applications.

## INTERPRETATION OF RESULTS

### For visual fluorescence detection

#### (For HumaLoop M)

Set each reaction tube into the Fluorescence Detection Unit, irradiate and observe the tube from the side.

#### (For other incubator using UV light)

Irradiate the bottom of each reaction tube and observe from the side through goggles or other UV-protective eye mask.

For a valid run, the following results must be obtained when read at the specified time:

- Positive control: Green fluorescent light is emitted.
- Negative control: No fluorescent light is emitted.

If any control is invalid, all samples in the run should be reported as invalid and the samples should be retested.

After confirming that the run is valid, evaluate samples as follows:

- Positive sample: Green fluorescent light is emitted.
- Negative sample: No fluorescent light is emitted.

#### For real-time turbidity detection with HumaTurb C+A

After confirming that the turbidity increases in the positive control but not in the negative control solution, evaluate samples in accordance with the following criteria (Fig 1 and 2).

- Positive: Some increase is observed in turbidity.
- Negative: No increase is observed in turbidity.

#### Notes:

- (1) The minimum detection sensitivities of the MALARIA Pf Detection Kit are 7.5 copies per test, respectively. In the case of a negative test result, patients with any persisting or worsening symptoms should be considered for retesting and other possible causes of symptoms should also be considered and investigated. The LAMP assay is highly sensitive and may detect low-level parasitaemia that are not the direct cause of the presenting symptoms. The patient's clinical condition should always be taken into account when making a final diagnosis and determining management.
- (2) Although the primers have been designed to target a region containing a relatively small number of variations, it is conceivable that *Plasmodium falciparum* infections may occur with variations in this region that are not well detected by this product. Therefore, a negative

test does not always rule out an infection by Pf.

- (3) This product is a kit for a qualitative detection; it is not designed for a quantitative measurement. The intensity of the fluorescent light observed or the rise time of turbidity measured by the Real-Time Turbidimeter does not correlate with the concentration of template DNA.

## INTERFERING SUBSTANCES

Our in-house studies have revealed that the turbidimetry measurement was not affected by the presence of EDTA (6.0 nmol/test), free bilirubin (1.2 µg/test), conjugated bilirubin (1.2 µg/test), chyle (formazine turbidity:14,400) hemolytic hemoglobin (3.4 µg/test) and human genomic DNA (1.2 µg/test). EDTA may cause false positive results when the result is read by fluorescence.

With regard to drugs, our in-house studies have revealed that the measurement was not affected by the presence of isoniazid (840 ng/test), ethambutol (204 ng/test), rifampicin (958 ng/test), pyrazinamide (4.2 µg/test), acetaminophen (322 ng/test), loxoprofen sodium (637 ng/test), cefotaxime sodium (12 µg/test), levofloxacin (270 ng/test), atovaquone (1.59 µg/test), proguanil (20 ng/test), chloroquine (40 ng/test), quinine (960 ng/test), doxycycline hydrochloride (360 ng/test), mefloquine (168 ng/test), primaquine (18 ng /test) and artemisinin (93 ng/test).

## PERFORMANCE CHARACTERISTICS

### (1) Sensitivity and accuracy

In testing the following samples:

- Negative sample (concentration: 0 copy/test)
- Positive sample 1 (100 copies/test)
- Positive sample 2 (1,000 copies/test);

The negative sample shall test negative and the positive samples 1 and 2 shall test positive.

### (2) Within-run reproducibility

In testing five negative and five positive samples simultaneously, the negative samples shall test negative throughout and the positive samples shall test positive throughout.

### (3) Limit of detection

7.5 copies/test

### (4) Cross-reactivity

The measurement system tested positive for *Plasmodium falciparum* and negative for other pathogens, as detailed in the table below:

<i>Plasmodium</i> genus	
<i>Plasmodium falciparum</i>	Positive
<i>Plasmodium vivax</i>	Negative
<i>Plasmodium ovale</i>	Negative
<i>Plasmodium malariae</i>	Negative
Other pathogens	
<i>Trypanosoma brucei</i>	Negative
<i>Trypanosoma cruzi</i>	Negative
<i>Leishmania donovani</i>	Negative
<i>Schistosoma mansoni</i>	Negative
<i>Theileria parva</i>	Negative
<i>Mycobacterium tuberculosis</i>	Negative
Influenza A virus	Negative
Human immunodeficiency virus	Negative

### (5) Information about a calibrator

The performance test for this product used plasmid DNA containing the mitochondrial DNA of Pf as a calibrator.

### (6) Clinical performance

Malaria is caused by certain parasites of the *Plasmodium* genus, which are transmitted via the bites of infected mosquitoes. After a period of development in the liver, blood stage parasites are released. The blood-stage parasites enter red blood cells, lysing the cells during subsequent reproduction and causing symptoms including fever. The LAMP assay detects DNA from blood-stage parasites.

The Malaria LAMP was evaluated in blood samples sent for malaria testing in a specialist laboratory. Following standard microscopy, samples were anonymised, DNA was extracted from the blood samples using the PURE<sup>6)</sup> method and tested by the Malaria LAMP kit using Pf-specific primers, with two different read-out methods. Finally, nested PCR was deployed as a reference method. 705 Malaria LAMP tests, of which 56 were from patients positive by microscopy, were compared to the nested PCR as a reference method. Sensitivity and specificity of the LAMP Malaria Pf were 98.4% and 98.1%, respectively.

	Sensitivity (95% CI)	Specificity (95% CI)	PPV*	NPV*
Malaria LAMP vs nested PCR	98.4%	98.1%	83.8%	99.8%
Malaria LAMP vs microscopy	100.0%	96.6%	70.3%	100.0%
	93.2 – 100.0	94.9 – 97.9		

\*positive and negative predictive values

## ORDERING INFORMATION

Product Code	Product Name	Contents
978000	<b>Loopamp™ MALARIA Pf Detection Kit</b>	<b>96 tests</b>
962000	HumaLoop M	1 Main Unit 1 Fluorescence Detection Unit
970000	Loopamp™ PURE DNA Extraction Kit	90 tests
963200	HumaTurb C+A	1 Control Unit 1 Amplification Unit
964000	HumaHeat	Heating block
980000	HuMax ITA	Bench-top centrifuge

## REFERENCES

- 1) Notomi T., et al.: Nucleic Acids Research 28, No. 12, e63 (2000)
- 2) Nagamine K., et al.: Clin. Chem. 47, No. 9, 1742–1743 (2001)
- 3) Mori Y., et al.: Biochem. Biophys. Res. Commun. 289, No. 1, 150–154 (2001)
- 4) Tomita N., et al.: Nat. Protoc. 3, No. 5, 877–882 (2008)
- 5) The guideline for the bio-safety and bio-hazard (by the Japanese Society for Bacteriology): Japanese Journal of Bacteriology 54, No. 3, 667–715 (1999)
- 6) <https://www.finddx.org/wp-content/uploads/2016/04/SOP-LAMP-Malaria-Aug2012.pdf>

## Flow chart

### Operation procedure for real-time turbidity detection

Preparation of sample solution

Prepare sample DNA solutions by extracting DNA from collected blood samples.

Reagent preparation

Take the required number of reaction tubes of dMAL Pf.

(For samples, negative and positive controls)

Mixing of reagent and sample solution

Transfer 30 µL of sample or control solution into each reaction tube.

Use NC Mal as negative control.  
Use PC Mal as positive control.  
(The positive control should be prepared at last.)

Invert the reaction tubes to collect the solution in the cap. Leave the tubes standing upside down for 2 minutes.

Invert the reaction tubes 5 times to mix the contents and then spin them down.

(Avoid making any bubbles.)


Amplification

Load the reaction tubes onto the reaction block of the turbidimeter.

As directed in the instructions for use of the turbidimeter, start reaction and measure and evaluate the turbidity (for 40 minutes at 65.0°C).

Confirm the completion of polymerase inactivation (for 5 minutes at 80°C, or for 2 minutes at 95°C). Take all reaction tubes from the turbidimeter and discard them without opening. Be careful not to damage the tubes.

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**EIKEN CHEMICAL CO., LTD.**