

Specification of OC-SENSOR iO	
Name	OC-SENSOR iO
Principle	Latex agglutination immuno-turbidimetry
Method	One step rate assay
Working Speed	88 tests / hour
Sample Rack	20 samples : 5 samples x 4 racks
Reaction Cuvette	Disposable acrylic cuvette 10 serial cell cuvette, Product Name 'DISPO-10'
Sampling System	Drawn up by sampling nozzle
Reagent Dispense	Drawn up by reagent nozzle
Mixing System	By spinning mixer
Thermostatic System	37°C Block Heater
Light Source	L.E.D. (wavelength 660nm)
Light Detector	Silicon Photodiode
Data Input	Color LCD touch panel (5.7 inches)
Memory Capacity	Sample: 5,000 tests, QC: 1,000 tests, Error: 1,000 tests
Data Output	Built-in thermal printer, RS-232C, USB
Barcodes	Sample barcode, Rack barcode
Power Required	AC100-240V 50/60Hz, 150VA
Dimension	W360 x D560 x H425 (mm) / 35kg

Reagents & Consumables



OC-Auto3 Latex Reagent / OC-Auto3 Buffer




OC-Control LV1 / OC-Control LV2
OC-Calibrator 2 kit



DISPO-10

SQUEEZE UP SYSTEM



The filtered sample goes into the sample cup part when squeezed by the instrument, then sampled into the reaction cuvette.

- Before use, please see attached operation manual.
- Above specifications may change in order to improve machine performance without notice.



MANUFACTURER


EIKEN CHEMICAL CO., LTD.
4-19-9, Taito, Taito-ku, Tokyo 110-8408, Japan
<http://www.eiken.co.jp>

3034 CP2

March, 2016

EIKEN CHEMICAL developed the first generation of "OC-SENSOR" in 1989.
 Integrating all the accomplishments for more than 20 years,
 We triumphantly introduce "OC-SENSOR io".

OPERATION PROCEDURE



Set the Sample racks.

Set the cuvettes (DISPO-10).

Set the Latex Reagent.

Press START to start analyzing.

The first test result will be printed out in 7 minutes.

OC-io	2010.08.02	Current date
*****		Calibration Curve No.
CC No.	1	SEQ No. (001-999)
001	001-01	Sample Rack No.
Y40MUU		Position No. in the rack
	28ng/mL	Flag (- / +)
002	001-02	Sample Barcode Number
T4KFEE		
	+ 345ng/mL	Concentration
003	001-03	Over Range (1001 over)
T7UESA		
	+OR2223ng/mL	Error Code
004	001-04	
	NO SMP L SOL	
005	001-05	
	- 0ng/mL	

- Compact Body**
W360 x D560 x H425 (mm) small size for small space.
- Simple Operation**
Just press "START" to start analyzing.
- Continuous Loading**
During analysis, continuous loading samples, latex reagent and cuvettes is available.
- High Throughput**
88tests / hour