Reagents

Faecal Immunochemical Test (FIT)







OC-FIT Control LV1 OC-FIT Control LV2 OC-FIT Control LV3 OC-FIT Calibrator

Faecal Calprotectin Test (FCa)



OC-FCa Reagent




Dimensions


OC-FCa Control LV1 OC-FCa Control LV2 OC-FCa Control LV3 OC-FCa Calibrator

All reagents are ready-to-use. OC-SENSOR Sample Diluent is also available.

Specifications

Product Name	OC-SENSOR PLEDIA	
Principle	Latex agglutination immunoturbidimetry	
Throughput	320 tests/hour	
Sample Loading Capacity	10-sample rack \times 20 *200-sample rack-tray (optional)	
STAT Sample	Up to 10 samples simultaneously	
Sample Dilution	15 times and 250 times	680
Reaction Cuvette	Semi-disposable acrylic cuvette (11 serial cells × 5 cuvettes) Automatic cuvette cleaning with wash solution and purified water	
Sampling System	Drawn up by sample nozzle (Liquid surface detection, automatic cleaning with wash solution and purified water)	
Reagent Dispense	Drawn up by reagent nozzle (Liquid surface detection, automatic cleaning with purified water)	
Mixing System	By spinning mixer (automatic cleaning with wash solution and purified water)	
Thermostatic System	Reaction wheel: air bath block heater, Reagent housing: block heater	
Light Source	LED (wavelength 660 nm)	
Light Detector	Photodiode	
Data Input	Color LCD touch panel (10.5 inches)	
Memory Capacity	100,000 test results	
Data Output	Built-in thermal printer, RS-232C, USB	
Barcodes	Rack barcode, sample barcode, reagent barcode, two-dimensional barcode (optional)	560
Dimensions	W630mm × D630mm × H560mm	
Weight	58kg	
Power Required	AC100-240V 50/60Hz, 500VA	
Applicable standards and laws	The PMDAct, CE, TGA, TFDA	






630 mm





Please use this product after reading the Operation Manual carefully.

Scan to learn more about OC-SENSOR PLEDIA



Tower Business Centre, 2nd Flr., Tower Street, Swatar, BKR 4013 Malta

CE IVD







FEATURES & PERFORMANCE



High Throughput Sampling

Improved processing speed, increased sampling loading capacity



Throughput

Capable of analysing 320 tests per hour.



Set reagents

With the number of reagent slots increased from two to three. PLEDIA is more efficient for high throughput population screening.



Onboard sample capacity

Up to 200 samples can be loaded at a time. Operation with a tray is also available. *The tray system is optional.

Sample rack detection sensor


Analysis starts automatically when the sample rack is set.



Advanced Usability Improved screen operability







User-friendly touch panel

10.5 inch LCD touch panel. The adjustable display angle improves visibility.

Real-time information

Real-time information is available, including information on reagents and sample processing time.

Quality control function

X-R control screen, CV auto calculation, etc. Massive QC function is available.

Random Access

Multiple Items Ordering

Bidirectional communication between PLEDIA and the host allows for multi-item order processing in each sample.









Fully Automated System

1 an 1 1 an 1 1 an 1 1 an 1 1 an 1 1 an 1		1 2 2 3 3 3 3			11111	
1		010	***		01	
1 an 1 7 an 1	0 444 0 1444	014	-		0	
1 an 1	-	010				
-			-		01	
	Cabe	-				
		- Cran	itek		0-	
3 ee 4	0-++	OPP	æ / .		01	
		11	1	11	10	1-
		<u> </u>	1.	lor	1	
	• •	• s s / •	e e e 7 7 e e	e e e 2 e / e = =		e e e e e e z e e e ee z e e e ee



Equipped with the capability for automatic startup, start, and one-touch shut down.

utomatic identification of easurement modes

Barcodes on racks enable the measurement mode to be identified automatically.

After completing an analysis


Retest racks enable simple re-testing.

With a dilute rack, 15-times and/or 250-times dilution (s) is automatically performed.



Dilute retest rack

Prozone recognition




If the \triangle Abs at the T6-8 photometric point of samples under the upper limit of measurement is above the lower limit of the PRC, the sample is recognized as within the prozone (PRC region).

Prozone of highly-concentrated samples is checked through the Primary Rate Check (PRC) method.

Squeeze-up system





As sampling bottle is squeezed, samples are filtered into the sample cup part, and then drawn by the sample nozzle.