## QUALITY CONTROL

Good laboratory practices recommend the use of appropriate controls to ensure the proper performance of the kit.

### Procedural Control

The Procedural Control is found in the Procedural Control Region of the test strip. This control assures the operator that (A) sample addition and migration through the test strip has occurred and that (B) the control anti-mouse antibody and the reporter MAb are intact and functional. This control does not ensure that the capture antibody is accurately detecting the presence or absence of hHb in the sample.

### External Control

External controls are used to assure the operator that the capture and conjugated antibodies are present and reactive. External controls will not detect an error in performing the patient test procedure. Controls should be assayed once per kit. To use, unscrew the white cap on the sample bottle. Add commercially available control with known hHb concentration. Replace the white cap and shake vigorously. Follow step three to five of the test procedure. If controls do not perform as expected, do not use the test results and repeat the test.

The buffer in sampling bottle may be used as a negative control.

## LIMITATIONS

- OC-Light™ S FIT is intended only for the detection of hHb in faeces. It is not advised for use in patients suspected of upper GI bleeding.
- Patients with the following conditions should not be considered for testing as these conditions may interfere with the test results:
  - Bleeding haemorrhoids
  - Constipation bleeding
  - Urinary bleeding
  - Menstrual bleeding
- Certain medications such as aspirin and non-steroidal anti-inflammatory drugs may cause gastrointestinal irritation and subsequent bleeding in some patients and cause positive results.
- As with any occult blood test, results obtained with OC-Light™ S FIT should not be considered conclusive evidence of the presence or absence of GI bleeding or pathology. The OC-Light™ S FIT is designed for preliminary screening. It is not intended to replace other diagnostic procedures such as colonoscopy or sigmoidoscopy in combination with double contrast barium x-ray.
- Because gastrointestinal lesions may bleed intermittently and blood in faeces is not distributed uniformly, a negative test result does not assure absence of lesion.
- Urine and excessive dilution of samples with water from the toilet bowl may cause erroneous test results.
- OC-Light™ S FIT is not for use in testing urine, gastric specimens or other body fluids.

## EXPECTED VALUES

Positive rates with immunochemical faecal occult blood tests have been shown to vary in each patient population depending on the test used, age and ethnicity of the patient and the predisposition to colorectal disease and other factors that may be associated with lower gastrointestinal bleeding. The OC-Light™ S FIT detects Hb in faeces at levels as low as 50 ng/mL or 10 µg/g faeces.

## PERFORMANCE CHARACTERISTICS

### Sensitivity

The sensitivity of OC-Light™ S FIT is 10 µg hHb/g of stool, or 50 ng hHb/mL of buffer.

The ability of OC-Light™ S FIT to detect hHb variants was determined by testing HbS and HbC, in comparison to a reference, HbA0. OC-Light™ S FIT detected the presence of all variants, HbA0, HbS and HbC.

### Interference Testing

Cross-reactivity studies were performed by adding the following non-human haemoglobin (Hb) and tissue extracts to OC-Light™ S FIT. Hb from bovine, equine, goat, porcine, rabbit, sheep, turkey, and fish was added to normal stool extracts containing 0 and 50 ng/mL hHb. All tests with 0 ng/mL hHb with non-human Hb were negative, indicating no cross-reactivity, and all with 50 ng/mL hHb were positive, indicating no interference to the test. The same test was performed with addition of tissue extracts from beef, chicken, fish, horse, pork, rabbit, goat, and sheep, resulting in no cross-reactivity or interference to the test.

### Dietary Testing

A potential interference of dietary substances on the OC-Light™ S FIT was assessed. Aqueous extracts of raw broccoli, cauliflower, cantaloupe, horseradish, red radish, parsnip and turnip were added to the test device to determine if vegetable extracts cross react with the test. The extracts were prepared by homogenizing raw vegetable in a food processor and then subsequently centrifuging the extract to separate the solid and liquid phases. Dietary Iron and Vitamin C supplements were also tested for cross-reactivity. No cross-reactivity was evident.

### Comparison Study

OC-Light™ S FIT was compared with a commercially available device, OC-LIGHT using 147 specimens. Human haemoglobin-free stool suspension was spiked with a known level of hHb to result in the following concentrations: 0 ng/mL, 25 ng/mL, 37.5 ng/mL, 50 ng/mL, 62.5 ng/mL, 75 ng/mL and 2000 ng/mL. The total number of determinations per level of hHb was 21. The overall agreement rate between OC-Light™ S FIT and OC-LIGHT was 99.3%, with positive agreement rate of 98.8% and negative agreement rate of 100%, demonstrating that the analytical performance of the device is substantially equivalent.

## PRODUCT CODE, PRODUCT NAME AND STORAGE

Product code	Product name	Contents	Storage
V-PC52	OC-Light™ S FIT Test Strips	10 x 50 strips	2-30 °C
V-PH82	OC-Light™ S FIT Sampling Bottles	20 x 50 bottles	2-30 °C

## REFERENCES

1. J.S.Mandel et al.: New England Journal of Medicine, 328: 1365-1371, 1993.
2. S.J.Winawer et al.: Journal of National Cancer Institute, 85 (16): 1311-1318, 1993.
3. N.Hiwatashi et al.: Japanese Journal of Cancer Research, 84: 1110-1112, 1993.
4. H.Saito et al.: International Journal of Cancer, 61: 465-469, 1995.
5. B.Levin et al.: Gastroenterology, 134(5): 1570-1595, 2008.
6. S.Ciatto et al.: British Journal of Cancer, 96: 218-221, 2007.
7. G.Grazzini et al.: European Journal of Cancer Prevention, 13: 16-26, 2006.
8. E.Fenocchi et al.: European Journal of Cancer Prevention, 15: 384-390, 2006.
9. T.Takesita et al.: Journal of Coloproctology, 38: 780-783, 1985.
10. LGM Rossum et al.: British Journal of Cancer, 101: 1274-1281, 2009.

Advena Ltd.

Tower Business Centre, 2nd Flr., Tower Street,  
Swatar, BKR 4013 Malta

EIKEN CHEMICAL CO., LTD.

4-19-9, Taito, Taito-ku, Tokyo 110-8408, Japan

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341760A-B

Español

## OC-Light™ S FIT Test Strips Sampling Bottles

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